
Neo-Safe® T CU 380 / CU 380 Mini

Intrauterinapparat zur Empfängnisverhütung
Intrauterine Contraceptive Device
Dispositivo Contraccettivo Intrauterino

DE

Neo-Safe® T CU 380 Neo-Safe® T CU 380 Mini

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PHYSICIAN INSTRUCTIONS


Presentation

The Neo-Safe® T CU 380 / Neo-Safe® T CU 380 Mini Intrauterine Device (IUD) is made of polyethylene and barium sulphate. Copper wire is wound around the vertical arm. The copper surface is 380 mm².

Dimension: The Neo-Safe® T CU 380 is 31,8 mm ± 0,2 mm wide and 31,9 mm ± 0,2 mm long.

The Neo-Safe® T CU 380 Mini is 24,0 mm ± 0,2 mm wide and 30,0 mm ± 0,2 mm long.

The polyethylene frame and the wire are radio-opaque.

- | | |
|------------------|-------------------------------|
| a) side arms | e) centimetre scale |
| b) vertical arm | f) insertion rod |
| c) nylon threads | g) marked part on the plunger |
| d) blue flange | h) plunger |

Characteristics

Neo-Safe® T CU 380 / Neo-Safe® T CU 380 Mini IUD offers almost complete protection against pregnancy. Correctly inserted, the IUD is safe for women at low risk of sexually transmitted disease.

The Neo-Safe® T CU 380 / Neo-Safe® T CU 380 Mini is an intrauterine contraceptive device made of polyethylene and barium sulphate. The banded side arms are flexible and shaped in such a way as to keep the IUD adjacent to the fundus, even in case of uterine contractions, without stretching the uterine cavity or touching the entrance of the fallopian tubes.

Neo-Safe® T CU 380 / Neo-Safe® T CU 380 Mini IUD prevents pregnancy by blocking fertilization.

The contraceptive effect is likely due to a pronounced sterile inflammatory reaction, which takes place as a result of a foreign body response in the uterus. The concentration of various types of white blood cells, prostaglandins and enzymes in uterine and tubal fluids increases markedly. These changes interfere with the movement of sperm in the genital tract, reducing their potential fertilizing capacity, so that fertilization is not possible. In the fallopian tubes, where fertilization is thought to take place, fewer sperm are found in copper IUD users than in non-users. Thus, the primary action is most likely altering the function or survival of sperm and ova before they can meet.

Indication

Female contraception (long lasting).

Contraindications (absolute)

1. Malignant diseases of the genital tract
2. Vaginal bleeding
3. Pregnancy
4. Past history of ectopic pregnancy or predisposing factors
5. Infections of the genital tract
6. Sexually transmitted diseases during the last 12 months (except bacterial vaginitis, repeated herpes infection, hepatitis B)
7. Abortion with infection during the last 3 months, pelvic inflammatory disease (PID)
8. Uterine malformations (congenital or acquired)
9. Allergy to copper

Contraindications (relative)

1. Valvular heart disease
2. Anaemia
3. Coagulation disorders
4. Anti-inflammatory treatment
5. Wilson's disease
6. Multiple sexual partners
7. Nulliparity

Warnings and precautions

Before inserting the IUD, a thorough medical history and an examination of the pelvic and abdominal cavity as well as a cervical smear are mandatory. Pregnancy, genital infections or sexually transmitted diseases should be excluded. The position of the uterus and the size of the uterine cavity must be determined to ensure correct insertion of the IUD. Neo-Safe® T CU 380 is designed for women with a uterine cavity depth of 6 - 9 cm. Neo-Safe® T CU 380 Mini is designed for women with a uterine cavity depth of less than 7 cm.

One month after the insertion of an IUD the woman must be re-examined to determine whether the IUD is properly placed and if there are signs of infection. Subsequent follow-up examinations are to be performed annually or more frequently if clinically indicated.

Pelvic inflammatory disease during IUD use should be treated without delay. For this reason the user must be instructed to report to her physician if there are suspicious signs such as vaginal discharge, pelvic pain and fever. This can be ascertained by gynaecological examination and/or ultrasound (if available).

If a treatment with antibiotics is not effective after 48 hours and signs of PID persist, the IUD must be removed immediately. Excessive bleeding or dysmenorrhoea during the first cycle after insertion should also be carefully assessed to ascertain if this is caused by the IUD, in which case it might have to be removed.

The possibility of perforation of the uterus during the insertion should always be considered, especially if the nylon threads are invisible or cannot be drawn out of the cervical canal. If there are any doubts about the position of the IUD (e.g. if the IUD is extremely difficult or painful to insert) the appropriate diagnostic techniques should be used (flat x-ray of the pelvis, ultrasonography, hysteroscopy, laparoscopy).

If the nylon threads appear to be longer than at insertion, an ultrasound should be carried out to determine if the IUD has been displaced, which might decrease its contraceptive efficacy.

If the threads are not visible at follow-up examination, pregnancy must be excluded. The threads may have been drawn up into the cervical canal or the uterus and usually reappear during the next menstruation. In case of doubt the IUD can be identified by feeling carefully with a suitable instrument, by ultrasound or by an x-ray after excluding pregnancy.

In case of accidental pregnancy with the IUD in situ, it is mandatory to determine (by ultrasound) whether the pregnancy is intrauterine or ectopic. Up to the end of the first trimester, if the threads are visible, the IUD should be carefully removed. After that the patient should be offered the option of elective abortion as soon as possible, since there is an increased risk of PID and many other obstetric problems, e.g. premature labour, placenta praevia and abruptio placentae.

However, if the woman wishes to continue her pregnancy, close monitoring is mandatory. She should be informed about the risks of keeping the IUD in situ. There are no contraindications to breastfeeding with an IUD in situ.

Perforation

Perforation or penetration of the uterine corpus or cervix by the IUD may occur, usually during insertion. The risk is increased during the postpartum period and slightly increased if the woman is fully breastfeeding. The device must be removed as soon as possible if this occurs.

Interactions

The available experience indicates that, in general, drug interference with the contraceptive efficacy of Neo-Safe® T CU 380 / Neo-Safe® T CU 380 Mini is highly unlikely. However, published reports appear to show diminished efficacy with long term use of non-steroidal anti-inflammatory drugs (especially acetylsalicylic acid) and of corticosteroids. Short term use in the treatment of dysmenorrhoea with non-steroidal anti-inflammatory drugs does not appear to reduce contraceptive efficacy. Do not perform diathermy (short wave and microwave) of the sacral or abdominal region since heating may cause heat injury to the surrounding tissue.

Side effects

During use of IUDs heavier menstrual bleeding, slight intermenstrual bleeding, anaemia, dysmenorrhoea, pain in the lower abdomen or sacral area as well as pelvic infection may occur. The risk of an ectopic pregnancy is increased. In rare cases, the IUD may penetrate into the wall of the uterus or perforate it. Allergic skin reactions may occur.

Physician instructions to IUD users

Users should be informed in detail as to the advantages and disadvantages of IUD contraception, not only so that they understand how it works, but above all so that any complications can be detected early. The user must learn how to feel the threads emerging from the cervical canal.

Insertion

First eliminate the possibility of pregnancy!

The best time for insertion is during menstruation to prevent insertion during non-diagnosed pregnancy. At this time the external and internal cervical os are physiologically dilated. This facilitates the insertion of the IUD without the need to dilate the canal in most instances. Alternatively, the IUD may be introduced within 5 days of unprotected coitus or within 15 minutes of delivery of the placenta or abortion in the first trimester. In these last two cases a higher expulsion rate is expected. Insertion immediately after unprotected coitus can increase the risk of PID. If the IUD cannot be inserted immediately after delivery of the placenta or abortion, insertion should be delayed for at least six weeks. In case of caesarean section insertion should be delayed for 12 weeks after delivery, although this is controversial.


Prior to insertion the vagina, cervix and cervical canal should be cleansed with an antiseptic solution. It is essential to determine the exact position of the uterus by bimanual pelvic palpation so that the Neo-Safe® T CU 380 / Neo-Safe® T CU 380 Mini can be inserted along its longitudinal axis. This can be accomplished by grasping the anterior or posterior lip of the cervix, depending on whether the uterus is anteverted or retroverted.

A local anaesthetic may be applied as a gel or can be injected in and around the cervix prior to insertion. Hysteroscopy should be carried out and the ascertained depth of the uterus be marked on the Neo-Safe® T CU 380 / Neo-Safe® T CU 380 Mini inserter tube with the aid of the movable flange.

Step-by-Step loading instructions


The IUD may be used by trained medical staff only. In order to minimize the risk of contamination we recommend using sterile gloves.

1. Open the sterile packaging of the Neo-Safe® T CU 380 / Neo-Safe® T CU 380 Mini halfway. Grasp the nylon threads and pull the device gently into the insertion tube (fig. 1) until the knobs at the end of the side arms cover the opening of the tube. The knobs should not be pulled into the tube.

Note: The ends of the side arms of the Neo-Safe® T CU 380 / Neo-Safe® T CU 380 Mini must not remain bended for over five minutes within the insertion tube; otherwise they may not bend back.


2. Steadying the blue flange with one hand, pull the insertion tube until the lower rim of the flange indicates the previously sounded length of the uterus on the scale (fig. 2).


Holding the threads straight in the tube, place the plunger into the insertion tube. This prevents that the threads will be disarranged by inserting the plunger.

Prior to insertion the tube can be bent to conform to the position of the uterus. The bending must be performed whilst the device remains in the sterile package after placing the plunger into the insertion tube.


3. Ensure that the blue flange indicates the direction in which the side arms of the IUD will open in the uterus. Remove the loaded insertion tube from the package and insert it through the cervical canal into the uterus until the blue flange touches the cervical os (fig. 3).


4. Hold the plunger tight and pull the insertion tube downwards until the lower rim of the tube reaches the marked part of the plunger. The horizontal arms of the IUD are now released in the uterus (fig. 4). Observe that the distance between the flange and the cervical os is now about 1,5 cm.

5. Holding the tube and the plunger together, gently push both as far as the blue flange again touches the cervical os (fig. 5). The IUD is now in the correct position.

6. Holding the plunger steady, pull the tube down to the ring of the plunger. The IUD is now entirely released in the uterus (fig. 6). To prevent pulling the device from the fundal position, first remove the plunger while keeping the insertion tube steady and only then remove the insertion tube.

Cut the threads leaving 2-3 cm visible outside the cervix.

Removal
Neo-Safe® T CU 380 / Neo-Safe® T CU 380 Mini should be replaced after 5 years.

Again, the most appropriate time to remove the IUD is during menstruation since both the internal and external cervical os are fully dilated. Grasp the threads of the IUD with a tenaculum and pull it along the longitudinal axis of the uterus. Try to insert the forceps at the entrance to the cervical canal in order to grasp the vertical arm of the IUD as soon as it has passed the internal cervical os. This avoids excessive tension on the threads, which could cause them to break. While doing this, distract the patient's attention by asking her to cough, and then remove the Neo-Safe® T CU 380 / Neo-Safe® T CU 380 Mini with one firm tug.

Insertion and removal of the IUD could cause slight pain and bleeding or may precipitate a vasovagal attack or an epileptic seizure. After removal, inspect the Neo-Safe® T CU 380 / Neo-Safe® T CU 380 Mini to see that none of its arms have been left in the uterine cavity.

In case of a lost IUD or lost parts of an IUD in the uterine cavity either hysteroscopy or ultrasound or x-rays should be used to determine its location: curettage may be advisable. In very rare cases of uterine perforation laparoscopy may be needed.

Packaging

1 x 1 sterile Neo-Safe® T CU 380 / Neo-Safe® T CU 380 Mini
10 x 1 sterile Neo-Safe® T CU 380 / Neo-Safe® T CU 380 Mini

Pharmaceutical information

Each IUD is sterilized with ethylene oxide and is intended for single use only. Do not resterilize. Do not use if package is damaged or open. Do not use after the expiry date printed on the sterile package. After removal the Neo-Safe® T CU 380 / Neo-Safe® T CU 380 Mini is to be disposed of in accordance with the guidelines in force for dealing with potentially infectious material.

Shelf-life: Four years.

Storage: Store in a dry place, below 35°C. Protect from direct sunlight and moisture.

Content: Copper, Polyethylene, Barium Sulphate, Polyamide 6.

Incompatibilities: No known incompatibilities.

Nature and contents of container: The device with accessories has been packed in a heat sealed sterilised pouch made of Tyvek + PET/PE.

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